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


Advanced Research Infrastructure for Materials and Nanotechnology in QST

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～ Sharing of Advanced Synchrotron Radiation Equipments and Experimental Data ～

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A Spoke Organization of Advanced Research Infrastructure for Materials and Nanotechnology in Japan

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National Institutes for Quantum Science and Technology (QST)

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National Institutes for Quantum Science and Technology (QST), manages two contract beamlines at SPring-8, a large synchrotron radiation facility, and QST's Synchrotron Radiation Research Center is developing advanced synchrotron techniques. QST has been designated by MEXT Advanced Research Infrastructure for Materials and Nanotechnology in Japan (ARIM) as a spoke organization in the special area "materials using quantum and electronic control to perform innovative functions". QST makes available the techniques described below for researchers from universities, companies, and public research organizations. We support analyses of crystal structures, local structures, electronic states, by using scattering, diffraction and absorption of synchrotron radiation.

散乱 Scattering


共鳴非弾性X線散乱
Resonant inelastic X-ray scattering

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
Please do not hesitate to contact us if there is any way we may be able to help with solutions to your R&D needs. Consultations with scientists or technicians of QST are available at any time, free of charge. Please refer to the QST-ARIM website <https://www.qst.go.jp/site/arim-english/> for contact details and the latest information.

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Only non-proprietary research proposals are accepted in the ARIM program. Proprietary research proposals can be accepted under QST's shared-use program.

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
Calls for proposals are held in May for experiments in the second half of the fiscal year, and in November for the first half of the following fiscal year. Urgent research proposals are welcome at any time. Only one experimental equipment can be used under a single research proposal. If multiple equipments are necessary, please make an individual research proposal for each equipment. Please contact the person in charge of the relevant equipment before submitting a proposal. There are two types of user supports: "Research Collaboration" and "Technical Assistance". Please contact us about the differences. ARIM proposals (non-proprietary proposals) are reviewed by the JAEA and QST joint proposal review committee.

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After conducting the experiments, please submit a user report by the specified deadline, including details of the experiment and its results. The user report will be published on the ARIM web in the next fiscal year.

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Data obtained through shared use of the equipment can be structured and stored in an easily usable format. QST encourages users to register your data as much as possible. Structured data that has been registered will only be accessible to the research team during the embargo period. After the embargo period is over, the structured data will be shared within the ARIM project. After that, it may be shared for public use beyond the ARIM project.

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SPring-8 users who have performed a non-proprietary experiment are required to publish their results in the either form of ① a refereed journal article (including refereed proceedings and dissertations), or ② a SPring-8/SACLA Research Report, refereed by JASRI, or ③ a technical journal article approved by JASRI. At least one publication must be registered in the SPring-8 Publications Database within three years from the end of the half-year research term during which the experiment was carried out.

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After all experiments of your research proposal were conducted, a handling charge and a usage fee are incurred. If you do not register any data higher usage fee is applied. For the latest information of them, please refer to the QST-ARIM website: https://www.qst.go.jp/site/arim-english/.

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
- A nuclear resonance spectrometer that generates ultra-high-brilliant 14.4 keV Mössbauer γ -rays from synchrotron radiation enables local analysis of magnetic interactions, electronic states, and spin configurations of functional substances and materials.
 - Focused Mössbauer γ -rays allow microscopic analysis of tiny parts of a sample and study of material states under ultrahigh pressure, whose sample space is limited to around 10 μ m.
 - Crystal site-selective magnetic structure analysis using a diffraction method and local magnetism survey of magnetic thin film using grazing incidence Mössbauer spectroscopy can be used.
 - Ni is also observable. Consultation is required for other elements.
- Typical application: Chemical analysis of earth-related materials under ultrahigh pressure, study of local magnetism on surface/interface of functional magnetic thin films etc.

共鳴非弾性X線散乱装置 *Resonant inelastic X-ray scattering spectroscopy station*

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
- This device enables X-ray spectroscopy experiments with an energy resolution of about 0.1 eV for both incident X-rays and scattered (emitted) X-rays.
 - In-situ observation of electronic states of catalytic materials is possible by high-energy resolution fluorescence detected X-ray absorption spectroscopy (HERFD-XAS) and X-ray emission spectroscopy (XES).
 - Using a resonant inelastic X-ray scattering (RIXS) technique, momentum-dependent electronic excitations can be measured in a four-circle geometry.
 - Sample temperatures of 10–800 K are available.
- Typical applications: Observation of charge, spin, and orbital excitations in transition-metal compounds, Operando analysis of electronic states in catalysts and electrodes.

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- This device enables real-time and in-situ observation using X-ray diffraction, making it ideal for studying growth processes of semiconductor heterostructures, multilayer films, quantum dots, and nanowires.
 - The growth of semiconductor nitrides such as GaN and InN is possible by using molecular beam epitaxy (MBE) with RF generated atomic nitrogens (RF-MBE).
 - This supports the creation of high-quality electronic materials.
- Typical application: Real-time analysis of growth processes for semiconductor quantum dots and semiconductor multilayer films.

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- ・ *This device allows observation of changes in material state and reaction progress under high pressure and high temperature using energy-dispersive X-ray diffraction with white X-rays.*
- ・ *This device can generate pressures up to 100,000 atmospheres and temperatures up to around 2000°C*
- ・ *This supports high-temperature and high-pressure synthesis of materials.*
- ・ *It can also generate a high-pressure and high-temperature hydrogen atmosphere of up to 100,000 atmospheres and around 1000°C, which is utilized for exploring new hydrogen storage materials through in-situ observation of metal hydrogenation reactions under high-temperature and high-pressure conditions.*

Typical application: In-situ observation of formation processes of metal hydrides under high pressure

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Diffraction for rapid-acquisition pair distribution function measurements

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- ・ *This device rapidly performs the measurements necessary for deriving the atomic pair distribution function (PDF), which allows evaluation of deviations from the average structure and local structure of crystals.*
- ・ *Using high-energy X-rays up to 70 keV, this device can perform X-ray total scattering measurements with a maximum Q of 27 \AA^{-1} , allowing derivation of PDF distance correlations up to about 100 Å.*
- ・ *This device allows low-temperature measurements using an open-flow cryogenic gas cooler, as well in-situ observation under hydrogen and nitrogen gas atmospheres of less than 1 MPa*
- ・ *Measurements can be performed in several seconds to several minutes per measurement point.*

Typical applications: Hydrogen absorbing materials, negative thermal expansion materials.

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Diffraction for coherent X-ray diffraction imaging

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
- ・ *"This device allows non-destructive 3D visualization of the internal structures (such as strain distribution, defects, voids, and domains) of nanocrystals (with grain sizes ranging from tens of nanometers to several micrometers) that are difficult to observe with electron microscopy, as well as their shapes.*
 - ・ *This method uses coherent X-ray diffraction imaging based on Bragg reflections in the hard X-ray region.*
 - ・ *This supports evaluations of nano-size materials.*
 - ・ *Please consult for the sample environment control. Sample temperature is controllable from a room temperature to around 1100 °C.*
- Typical applications: Three-dimensional visualization of 100–500 nm barium titanate nanocrystals. Three-dimensional visualization of 40 nm palladium nanocrystals.

SPring-8へのアクセス Access to SPring-8

http://www.spring8.or.jp/ja/about_us/access/

http://www.spring8.or.jp/en/about_us/access/

JR線とバスでのアクセス Access by JR-WEST railroad lines and bus routes

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"Advanced Research Infrastructure for Materials and Nanotechnology in Japan" is a program sponsored by the Ministry of Education, Culture, Sports, Science and Technology (MEXT) of Japan.